

# BEYOND CATALOGUE



## BEYOND YOUR IMAGINATION



Bio-photonic Biocommunication & Bio-cybernetic - Bio-feedback Systems  
Proton & Electron Resonance Technology

“ BEYOND is a photon based light oscillator with advance bio-harmonizing functionalities offering superior signaling for re-synchronization of bio-regulatory processes

Beyond is a product of an advance technological know-how with extraordinary capacities and World's first quantum frequency emitter including 4 unique features in combined form of fractal frequency delivery capability.


Device has the capacity to make the most minute calculations through its unique measurement technology making required frequency measurements with advance mathematical precision delivering personalized application modality via self-adapting capability. Through its autofeedback and automatic real time testing function, Beyond, provides ideal optimized testing and application for users.

Bio - Stimulation

Bio - Feedback

Bio - Photonic


Bio - Regulatory



## 4 in 1

Beyond is a product of advance German technology and developed to provide most accurate application capacities considering both physical and mental/emotional factors. Beyond is the first bio-regulation equipment utilizing **4 different technologies in the World.**

Beyond is a friend, providing assistance and help in resolving **most distinct health issues**, literally, at the touch of your fingertips.



## UNIFIED LIGHT APPLICATION UNIT

BEYOND, has the capacity of recalculating the signals received from the medium with careful mathematical accuracy and comes up with precise coherent feedback for **reharmonization** increasing the regeneration capacity of the body. Thus, improving the **self-regeneration** processes by providing required frequencies based on Proton and Electron resonance.

## QUANTUM LIGHT & MAGNETIC RESONANCE

Beyond offers personalized magnetic quantum light application possibilities for users.



Re-enforced Aluminum Body for better durability – Where the German Technology & Esthetics Meet

HIGH QUALITY  
MADE IN GERMANY

## 4 DIFFERENT TECHNOLOGIES IN THE WORLD

Beyond is the first bio-regulation equipment utilizing 4 different technologies in the World.



BEYOND harmonizes aimed metabolic processes with utmost precision and accuracy through Red & Infra-Red Light, Bio-resonance, Cymatic Resonance and Crystal Therapy providing a profound health application using Proton & Electron Resonance Technology. It's world's first Quantum Light Oscillator calibrated based on Proton & Electron Resonance.

Through its extensive information exchange technology, BEYOND provides required harmonizing signals to sustain and support ongoing communication on cellular level increasing the regeneration capability of our body.

BEYOND, increases the ATP (energy) production of our cells effecting the intracellular processes such as ATP production of the Mitochondria improving the autoregulation and regenerative processes of our bodies.









## 71 DIFFERENT APPLICATION OPTIONS BE YOUR OWN DOCTOR!

Through its touchscreen feature and easy-to-use operation capability, Beyond, provides an easy access to desired programs for application among 71 different program options.

Beyond offers detailed explanations through integrated information panel for each and every program through its advance program guide for users to have required information in utilizing the device, based on personal need and request.



## CATEGORY SPECIFIC BUTTONS



Each category is highlighted in a different color and icon for better navigation

| BONE     | Bone Regeneration<br>Osteoporosis<br>Arthritis                       | Inflammation<br>Infections<br>Bacteriostatic<br>Intolerance | Pressure Headache<br>Forehead Relaxation<br>Falling Asleep |
|----------|--|---|--|
| MUSCLE   | Muscle Stimulation<br>Muscle Regeneration<br>Muscle Relaxation       | Allergy<br>Immunity<br>Vitalization                         | Deep Sleep<br>Narcolepsy<br>Learning                       |
| INTERNAL | Reduce Muscle Tension<br>Increase Muscle Tone<br>Muscle Coordination | Anti-fungal<br>Adaptation<br>Viral Infection                | Short Memory<br>Long Memory<br>Perception                  |
| IMMUNE   | Muscle Cramps<br>Neck Relaxation<br>Breathing                        | Cell Regeneration<br>Cell Metabolism<br>Oxygen              | Meditation<br>Stress Reduction<br>Neural Activity          |
| CELL     | Circulation<br>Blood Pressure<br>Heart                               | Blood<br>E-smog<br>Detoxification                           | Stuttering<br>Dysphonia<br>Auditory Disorder               |
| PAIN     | Spleen<br>Stomach<br>Kidney  | Schumann 1<br>Schumann 2<br>Schumann 3                      | Hearing Loss<br>Tinnitus<br>Vision                         |
| SLEEP    | Liver<br>Reproduction<br>Pancreas/Sugar                              | Schumann 4<br>Acute Pain<br>Chronic Pain                    | Sense of Touch<br>Serotonin<br>Endorphin CNS               |
| LEARN    | Ovaries/Prostate<br>Estrogen/Testosterone<br>Wound Healing           | Endorphin PNS<br>Tension Headache                           | Addiction<br>Neuropathy                                    |
| NERVES   |  |   |  |


## RUBY CRYSTAL - JEWEL OF THE LIGHT

Beyond offers a customized health application via  
**Ruby Crystal Emission Technology**

Enriched with **%100 Natural** Ruby Crystal

BEYOND also offers **Crystal Therapy** through Crystal  
Information Technology


91pcs Red & Infrared LED




## PERSONAL MOBILE QUANTUM LIGHT DEVICE!

BEYOND is World's first Mobile Light Equipment based  
on Proton & Electron Resonance Technology

Beyond technology allows better **Intracellular** and  
**Extracellular** communication supporting the cellular  
communication supporting the self-healing mechanism  
of the body.



Better Application with Atomic Precision



## SAFE AND RELIABLE WITH NO SIDE EFFECTS!

### Non-invasive



Due its low-frequency emissions, BEYOND has no side-effects and can be used safely. It does not produce nor provide any aggressive or invasive health applications.



### LONG LASTING BATTERY LIFE..

With its long-lasting battery life and fast charging capability BEYOND provides desired features for better utilization.



- ON / OFF Switch
- USB Socket
- Charging Socket
- Light Indicator





## 3 DIFFERENT COLOR OPTIONS THE COLORS THAT FIT YOUR TASTE



Metallic Silver



Champagne Gold



Metallic Rose

With its 3 different color options, BEYOND offers you the opportunity to choose the color option you desire according to your taste and offers 7 different language options.

## 7 DIFFERENT LANGUAGE OPTIONS

English  
Italian

French  
Turkish

German  
Russian

Dutch

## AESTHETIC DESIGN

Aesthetic design, elegant and minimal appearance with customized colors.



## AT THE TOUCH OF YOUR FINGERTIPS..


Beyond provides easy and fast access to all programs with its touchscreen. Easy software operation guide with **user friendly interface**.

BEYOND is equipped with state-of-the-art new generation integrated electronic components. The touchscreen allows users to self-manage and articulate the desired programs through its **4.3" TFT Display** with **480x240** resolution.



### ALUMINUM CASE

Durability and Longevity  
Reinforced Aluminum Body for better durability




## SAFE & ERGONOMIC

With its easy-to-use interface, BEYOND has designed in a way that everyone can use with no complication.


BEYOND is a device developed to be used everyone safely. It is a **"Wellness"** device suitable for both professional and home use.



## HOW TO USE BEYOND?




Please Switch-ON the device by pressing the ON/OFF button at the right-bottom of the device. The startup screen will appear within a few seconds after the switch ON and the device is ready for use.



Please place the device on the body where the face with LED Light Panel is in full contact with the skin. Place the device where the application will be made. For instance, your leg, your head, or any part of your body where the application is needed.



Please select the desired application program from the program list and press the PLAY button (Start) to initiate the initial testing.



After the 10-second initial testing and measurement, the device will automatically start the application.



10 consecutive beep sound will indicate that the application is completed.

## WATERPROOF PROTECTIVE COVER BAG

The protective cover bag of BEYOND is made of exclusive **premium quality fabrics** to ensure safety and required protection for the device. Cushiony interior surface of the cover bag provides **ultimate protection**.



The fabric the cover bag is made of provides **high durability** with scratch and **waterproof** capacities allowing users to carry their device anywhere they desire with absolute safety and protection.





## ELEGANT BOX DESIGN

---

The Box of BEYOND is carefully designed with best quality materials available to provide comfort and elegance



---

BEYOND is portable unit with extensive extraordinary features assisting your health, physically, mentally and emotionally.



beyond

L U X

Quantum Light Oscillator



ALC  
Cybernetic Technologies